


SPECIMEN FLOWCHART

VN-BECS V2.0 SYSTEM INTEGRATION

DOCUMENT: SPECIMEN PIPELINE
STATUS: **APPROVED**
REF: SOP 02 / SOP 08 / SOP 11



STEP	SUBSYSTEM / FEATURE	WORKFLOW & VERIFICATION LOGIC
01	Donor Center (LIMS)	PHLEBOTOMY & SAMPLE EXTRACTION When phlebotomy completes, bind the whole blood bag and sample tube using the ISBT-128 standard DIN bar code. Pre-dispatch is automatically registered and pushed to the core testing lab.

STEP	SUBSYSTEM / FEATURE	WORKFLOW & VERIFICATION LOGIC
02	Immunohematology & Reference Lab (IRL)	SPECIMEN RECEPTION & WORKLIST DISPATCH Specimen arrives at the lab. The control tower lists the sample as 'Pending Reception' (待接收檢體). Signing off automatically dispatches the tube into parallel workflows: serological IDM/NAT screening and crossmatch queues.
03	Analyzer Integration (LIS)	AUTOMATED TESTING & LIS SYNC The specimen runs inside automated analytical platforms (e.g. cobas 6800, ARCHITECT). High-throughput assay results for HIV, HBV, HCV, NAT, and Serology are synced directly back to the BECS database.
04	Quality Assurance (QA) Node	FINAL QUALITY RELEASE & INVENTORY UNLOCK Once all IDM, NAT, and blood grouping checks are reactive-free and fully compatible, the system triggers the 'Quality Released' status. This immediately unlocks the corresponding physical blood bags inside the Blood Supply Hub from TESTING to AVAILABLE status.